



Easy Pure RO (Reverse Osmosis)

CONTACT US

Phone: 310-830-6601

FAX: 310-830-9336

Email: info@sepor.com
& services@sepor.com

LOCATIONS

WEST COAST

718 N. Fries Avenue
Wilmington, CA 90744

EAST COAST

3740 NW 124th Ave
Coral Springs, FL 33065



Description:

The Easy Pure RO system is designed for pretreatment of lab water for deionization systems or purifying water for general laboratory use. To achieve Type I reagent grade water (Ultrapure), the water from the RO system must be further treated in a Deionization System. The Easy Pure RO is compact in size, has flow rates of 10 liters/hour, is microprocessor controlled, LED indicators, built in 6.5-liter product reservoir and has 24 hour standby. Easy Pure units have a prefilter cartridge and two membrane cartridges. Dimensions: 31 x 48 x 46 cm

Catalog Number	Description
130E-100	120v/1Ph/60Hz
130E-110	240v/1Ph/50Hz
130E-120	Easy Pure Cartridges